

1

2

3

4

5

6

0

Quote dwg

Date 3/30/09

A


B

C

A

B

C



Notes:

- 1. Tubing USP Class VI, Shore 60A
- 2. Packaging: Bulk packaged.
- 3. No Particulates or debris

Quote Drawing

Dimensional Tolerances
 All tubing dimension as specified.
 Fractional dimensions (+/- 1/16")
 0.0 (+/-0.1"), 0.09 (+/-0.05), 0.090(+/-0.005)

Optimum Processing Inc.
 336 Bon Air Center #128
 Greenbrae, CA 94904
 ph: 415.722.5747

21-1/4 x 4 Tubing Assy's

SIZE	SCALE	Third angle projection
A	NTS	

Drw: DM Smith **Date 3/30/09**
 This specification contains proprietary information. No reproduction of this document is allowed without approval or consent of Optimum Processing Inc.

21-1/4" Neoprene x 4" Silicone
 with Clamps, Plugs Tubing
 Assembly

Rev
 0

1

2

3

4

5

6