

1

2

3

4

5

6

0

Quote dwg

Date 3/30/09

A

A

B


B

C

C

D

D



Notes:

- 1. Tubing USP Class VI, Shore 60A
- 2. Packaging: Bulk packaged.
- 3. No Particulates or debris

# Quote Drawing

**Dimensional Tolerances**  
 All tubing dimension as specified.  
 Fractional dimensions (+/- 1/16")  
 0.0 (+/-0.1"), 0.09 (+/-0.05), 0.090(+/-0.005)

**Optimum Processing Inc.**  
 336 Bon Air Center #128  
 Greenbrae, CA 94904  
 ph: 415.722.5747

25mm x 5 in Filter Tubing Assy.

<b>SIZE</b>	<b>SCALE</b>	<b>Third angle projection</b>
A	NTS	

<b>Drw: DM Smith</b>	<b>Date 3/30/09</b>
----------------------	---------------------

This specification contains proprietary information. No reproduction of this document is allowed without approval or consent of Optimum Processing Inc.

1/8"x1/4"x5-1/2"L Silicone,.2um  
 25mm Filter Tubing Assembly

Rev

0

1

2

3

4

5

6