



# Quote Drawing

**Note:**  
1. Individually Single Bagged & Sterilized

**Dimensional Tolerances**  
All tubing dimension as specified.  
Fractional dimensions (+/- 1/16")  
0.0 (+/-0.1"), 0.00 (+/-0.05), 0.000(+/-0.005)

**Optimum Processing Inc.**  
336 Bon Air Center #128  
Greenbrae, CA 94904  
ph: 415.722.5747

1000ml EMF Assy OP

<b>SIZE</b>	<b>SCALE</b>	Third angle projection
<b>A</b>	<b>NTS</b>	

Drw: DM Smith Date 6/4/09  
This specification contains proprietary information. No reproduction of this document is allowed without approval or consent of Optimum Processing Inc.

**1L Corning EMF FlexVent Assy., Swabable Valve, Slide Clamp, .2u Filter**

Rev  
**1**